

VIVA

RF System






VIVA RF System

Auto Mode

- Automatic step up power output
- Initial power output of 5~200 watts, it increases by 5 watts every minute (0.08 w/sec) till the first roll off, followed by pulsing RF controlled by impedance based algorithm.



Continuance Mode

- Continuous power output
- The RF is generated continuously. Ideal mode for Thyroid, Osteoid Osteoma and track ablation.



Temperature Mode

- Continuous maintenance for target of temperature
- Automatically adjust the RF output to maintain and reach the target temperature



star RF Electrode

star RF Electrode-Fixed	Basic electrode
star RF Electrode-Injectable	Injectable
star RF Electrode-Uni_faced	Half-moon shape

VIVA RF Electrode

Adjustable



Accessories




Peristaltic pump




Single foot switch



Double foot switch



No.	Feature name	No.	Feature name
1	CONTINUANCE Mode button	5	RF START/STOP button
2	AUTO Mode button	6	Digital Screen
3	OHM/INFO button	7	GROUND PADS connectors
4	RF POWER control dial	8	RF CABLE connector



No.	Feature name	No.	Feature name
1	Ventilation fan	5	Foot switch connector (lower part)
2	Foot switch connector (upper part)	6	Equipotential grounding terminal
3	Data port	7	Fuse/power socket
4	Power switch		

Specifications

Rated power

Voltage	110/220 V~
Voltage range	100-230 V~
Maximum input voltage	250 V~
Maximum input power	450VA
Maximum input current	4A (100/120 V~ units), 2A (220/240 V~ units)
Fuse capacity	F5 AH 250V
Power frequency	50/60 Hz

Impedance measurement

Range	10-800 ohms
Resolution	1 ohm
Accuracy	10-50 ohms ± 10 ohm 51-300 ohms $\pm 15\%$ 301-800 ohms $\pm 30\%$

Radiofrequency output

Watts	0-200watts max output@ 50ohm
Accuracy	$\pm 20\%$
Resolution	1 watt
Frequency	480 kHz $\pm 10\%$
Drive on time	30minutes max.

Temperature measurement

Range	5°C-95°C
Resolution	1°C
Accuracy	± 5 °C

Operating Environment

Clean, dry area Temperature	15°C-40°C
Humidity	15-80% relative, non-condensing
Atmospheric Pressure	800-1060 hPa



www.STARmed4u.com

E-mail. inquiry@STARmed4u.com

Facebook. <https://www.facebook.com/STARmed4u>

ST-CA-12 | Rev.3 | May, 2019