











Electrodes with 3 Needles

- Effective ablation in close cluster of 0.5cm distance
- Easily separated as 3 individual electrodes for easy placement
- cost effectively eliminate the need to stock cluster type Electrodes and the 2 or 3 needles set
- Faster and Larger Ablation when used with VIVA Multi RF Generator System
- Positive user feedback on sharper Needle Tip for easier percutaneous approach

Octopus RF Electrode

Gauge	Length	Active tip
17G	15cm, 20cm	20mm, 25mm, 30mm, 40mm, 50mm
15G	15cm, 20cm	20mm, 25mm, 30mm, 40mm, 50mm

Octopus RF Electrode-Duo

Gauge	Length	Active tip
17G	15cm, 20cm	20mm, 25mm, 30mm, 40mm
15G	15cm, 20cm	20mm, 25mm, 30mm, 40mm

^{***}Optional selection of Octopus two electrode or three electrode

