

www.apro-korea.co.kr

CE 1639

## AUTOMATIC STOP

AK-B250  
B300  
B400



Electro-surgical Unit

# SENS

Lorem ipsum

Lorem ipsum

AK-B250  
B300  
B400



Ergonomic design and various functions

- Cut: pure cut / blend1 / blend2
- Coagulation: contact / spray
- Bi-polar: Imp. On & Off

Precise cutting, easy hemostasis, and no carbonized coagulation

Various clinical applications: General surgery, Gynecology, Dermatology, ENT, etc.

Safe and easy operation

- Automatic power stop before carbonization
- Patient plate monitoring for proper placement

Each function has each sound


**AK**  
APRO KOREA Inc.



MEDVAL  
healthcare



## SENS standard Accessories



1. Mono-polar Foot Switch
2. Bipolar Foot Switch
3. Disposable Plate Cable
4. Hand Piece
5. Standard Electrodes
6. Disposable Plate
7. Button Handle
8. Bipolar Cable
9. Metal Patient Plate

				AK-B250	AK-B300	AK-B400
<b>Input Voltage &amp; Frequency</b>				<b>AC 220V , 50/60Hz</b>		
<b>Power consumption</b>				<b>900VA</b>		
<b>Radio Frequency</b>				<b>Mono/400, Bipolar/500KHz</b>		
<b>Modulation Frequency</b>				<b>300KHz</b>		
<b>RF output power</b>	<b>Cut</b>	under 500 Ω	<b>Cut(pure)</b>	250W	300W	400W
		under 300 Ω	<b>Blend 1</b>	150W	230W	230W
		under 300 Ω	<b>Blend 2</b>	120W	180W	180W
	<b>Coagulation</b>	under 300 Ω	<b>Contact</b>	100W	120W	120W
		under 500 Ω	<b>Spray</b>	80W	100W	100W
	<b>Bi-polar</b>	under 100 Ω	<b>Imp. On</b>	80W	100W	100W
<b>Imp. Off</b>			80W	100W	100W	
<b>Power Step</b>				<b>1W</b>		
<b>Classification</b>				<b>CF(Cardiac Floating)</b>		
<b>Dimension(WxDxH)</b>				<b>341x394x160 mm</b>		
<b>Weight</b>				<b>16Kg</b>		